



# Product catalogue



---

# Contents

Protective apparel	4
Bed protection	20
Sampling & examining	24
Patient hygiene & care	30
Wound care	35
Surgical tapes	41
Adhesive tapes	42

## Protective apparel

Correct use of protective apparel prevents microbes and infection from spreading from personnel and visitors to patients and vice versa. Disposable PE-apron protects better than textile protective gown, where the particles are stored in the structure of the material, and is recommended to be used when there is a risk of working clothes getting wet or soiled. PE gloves give extra protection for example when used together with washing/examination gloves.

SELEFA® vinyl, latex and nitrile gloves offer a cost efficient solution both in health care and non-medical settings

Protective apparel assortment contains shoe covers, arm covers, caps and overalls etc. Many of these products have an important role in infection control but are also very well suited to non-medical uses.

SELEFA® assortment contains both medical devices (MDD 93/42/EEC) and personal protective equipment (PPE 89/686/EEC). All medical devices in this catalogue are Class I or Class Is (sterile) and all personal protective equipment are Cat I.

The labeling follows EN 980 and EN1041 standards. All products are equipped with barcodes to help handling in the warehouse.



Comfortably in control with  
SELEFA®

## Comfortably in control with SELEFA® aprons

Disposable aprons made of fluid resistant polyethylene are an important part of infection control strategy. They also help protect clothing from spills and splashes during medical procedures and cleaning tasks.

All Selefa aprons are in conformity with the European Directive for Personal Protective Equipment 89/686/EEC, Cat I.

Lightweight and disposable, these aprons are also ideal for use in food preparation and in industrial and pharmaceutical applications.

Selefa polyethylene aprons are available in a variety of sizes (widths/lengths), thicknesses and colours.

The packing systems available are

- Flat packed, 210405S
- On a block hanger, 210431
- On a roll in dispenser eg 210433
- Individually folded in a dispenser eg 21043050 (Bosco)
- Individually wrapped in a dispenser eg 210424

Available in SELEFA® assortment are handy apron holders made of plexiglass.



## Aprons without sleeves

### Apron

#### REF 210415

- White, LDPE film, embossed
- Size: 72 x 118 cm
- Packed in a handy dispenser box
- Pcs/box: 100
- Pcs/carton: 500



### Apron Bosco

#### REF 21043050

- Transparent, LDPE film, embossed
- Size: 81 x 125 cm
- Individually folded
- Packed in a handy dispenser box
- Pcs/box: 50
- Pcs/carton: 500



### Apron on a roll

#### REF 210433

- Transparent, LDPE film, embossed
- Size: 81 x 125 cm
- Packed into a handy dispenser box
- Pcs/box: 50 on a roll
- Pcs/carton: 500



### Apron on a block hanger

#### REF 210432

- Transparent, LDPE film, embossed
- Size: 81 x 125 cm
- Easy to dispense from hanger
- Pcs/bag: 50 on block hanger
- Pcs/carton: 500



### Apron

#### REF 210405S

- White, LDPE film, embossed
- Size: 81 x 125 cm
- Packed into a plastic bag with cardboard support
- Pcs/bag: 100
- Pcs/carton: 500



### Apron

#### REF 210406S

- White, LDPE film, embossed
- Size: 80 x 140 cm
- Pcs/bag: 100
- Pcs/carton: 500



## Long aprons without sleeves

### Apron, extra long

#### REF 210434

- Transparent, LDPE film, embossed
- Size: 81 x 150 cm
- Individually folded
- Packed into a handy dispenser box
- Pcs/box: 40
- Pcs/carton: 400



### Apron, extra long, in a roll

#### REF 210435

- Transparent, LDPE film, embossed
- Size: 81 x 150 cm
- Packed into a handy dispenser box
- Pcs/box: 50
- Pcs/carton: 500



### Apron, extra long

#### REF 210408S

- White, LDPE film, embossed
- Size: 85 x 160 cm
- Pcs/bag: 100
- Pcs/carton: 500



### Apron, extra long

#### REF 210407S

- White, LDPE film, embossed
- Size: 90 x 160 cm
- Pcs/bag: 100
- Pcs/carton: 500



## Heavy duty aprons without sleeves

### Apron, extra strong

#### REF 211400

- White, LDPE film
- Size: 85 x 140 cm
- Pcs/bag: 100
- Pcs/carton: 500



### Apron, extra strong, on a block hanger

#### REF 210431

- Transparent, LDPE film, embossed
- Size: 100 x 150 cm
- Extra wide
- Easy to dispense from block hanger
- Pcs/bag: 50 on block hanger
- Pcs/carton: 500



## Aprons with sleeves

### Apron with short sleeves

#### REF 210426

- Transparent, LDPE film, embossed
- Size: 98 x 120 cm
- Individually wrapped
- Packed into a handy dispenser box
- Pcs/box: 30
- Pcs/carton: 150



### Apron with sleeves, thumb loop

#### REF 210425T

- Transparent, LDPE film, embossed
- Size: 98 x 120 cm
- Individually wrapped
- Packed into a handy dispenser box
- Pcs/box: 20
- Pcs/carton: 100



### Apron with sleeves, thumb loop

#### REF 210425L

- Blue, LDPE film, embossed
- Size: 98 x 120 cm
- Individually wrapped
- Packed into a handy dispenser box
- Pcs/box: 20
- Pcs/carton: 100



### Apron with sleeves, thumb loop

#### REF 210424

- Blue, LDPE film, embossed
- Individually wrapped
- With thumb loop
- Size: 94 x 140 cm
- Individually wrapped
- Pcs/unit: 1
- Pcs/carton: 100






# SELEFA® PE gloves

## PE gloves, Extra Strong

- Transparent, smooth
- Available in 3 sizes
- Length 34 cm
- Pcs/bag: 100
- Pcs/carton: 5000

**REF 210329S**  
**REF 210330M**  
**REF 210331L**

Small  
Medium  
Large




Smooth surface

## PE gloves

- Transparent, smooth
- Available in 3 sizes
- Length 35 cm
- Pcs/bag: 100
- Pcs/carton: 5000

**REF 210320**  
**REF 210321**  
**REF 210322**

Small  
Medium  
Large




Smooth surface

## PE gloves

- Transparent, embossed
- Available in 3 sizes
- Length 33 cm
- Pcs/bag: 100
- Pcs/carton: 5000

**REF 21030333**  
**REF 21030533**  
**REF 21030433**

Small  
Medium  
Large




Embossed surface

## PE gloves

- Transparent, embossed
- Available in 4 sizes
- Length 34 cm
- Pcs/bag: 100
- Pcs/carton: 5000

**REF 210298**  
**REF 210299S**  
**REF 210300M**  
**REF 210301L**

X-small  
Small  
Medium  
Large




Embossed surface

# SELEFA® PE gloves

## PE gloves

- Transparent, embossed
- Available in 3 sizes
- Length 34 cm
- Packed in handy dispenser box
- Pcs/box: 400
- Pcs/carton: 4000

**REF 210299A** Small  
**REF 210300A** Medium  
**REF 210301A** Large




Embossed surface

## PE gloves, Extra Long

- Transparent, embossed
- One-size
- Length 53 cm
- Pcs/box: 100
- Pcs/carton: 2500

**REF 210312** Medium




X-long, embossed surface

## PE gloves, Extra Strong

- Rose, embossed
- Available in 3 sizes
- Length 35 cm
- Pcs/bag: 100
- Pcs/carton: 5000

**REF 10204** Small  
**REF 10206** Medium  
**REF 10205** Large




Embossed surface

## PE gloves, Extra Strong

- Rose, embossed
- Available in 3 sizes
- Length 34 cm
- Pcs/bag: 100
- Pcs/carton: 5000

**REF 210332S** Small  
**REF 210333M** Medium  
**REF 210334L** Large



Smooth surface

# Powder-free paper backed copolymer gloves

## Powder-free paper backed copolymer gloves, unsterile

- Available in 3 sizes
- Length 28 cm
- Transparent
- Smooth surface
- Pcs/dispenser: 100
- Pcs/carton: 1000

**REF 210319S** Small  
**REF 210319M** Medium  
**REF 210319L** Large



## Powder-free paper backed copolymer gloves, sterile

- Available in 3 sizes
- Length 28 cm
- Transparent
- Smooth surface
- Gamma-sterilized
- Pcs/unit pack: 1
- Pcs/dispenser: 100
- Pcs/carton: 2000

**REF 210318S** Small  
**REF 210318M** Medium  
**REF 210318L** Large



# SELEFA® cotton gloves

## SELEFA® Extra

- Natural white, strong
- Available in 3 sizes
- Length 30 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222884L S  
REF 222885L M  
REF 222886L L



## SELEFA® Special

- White
- Available in 3 sizes
- Length 30 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222889L S  
REF 222890L M  
REF 222891L L



## SELEFA® Standard

- White
- Available in 2 sizes
- Length 25 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222890 M  
REF 222891 L

## SELEFA® Consumer

- Natural white
- Available in 3 sizes
- Length 30 cm
- Pcs/box: 10
- Pcs/carton: 1600

REF 222895ADAS S  
REF 222895ADAM M  
REF 222895ADAL L



## SELEFA® Basic

- Natural white
- Available in 3 sizes
- Length 30 cm
- Pcs/bag: 40
- Pcs/carton: 2000

REF 222895S S  
REF 222895M M  
REF 222895L L



## Knitted cotton gloves with grip points

- Available in white and black
- Available in 4 sizes
- Length 22 cm
- Pairs/bag: 12
- Pairs/carton: 600  
(50 bags of 12 pairs)

White:  
REF 222901 S  
REF 222902 M  
REF 222903 L  
REF 222904 XL

Black:  
REF 222911 S  
REF 222912 M  
REF 222913 L  
REF 222914 XL



# SELEFA® vinyl gloves

## Vinyl

Powder-free SELEFA® vinyl glove offers practical protection for many applications. The 100 % synthetic glove is somewhat softer than the standard vinyl and has a good feel and performance. It is a good alternative for those suffering from Type I allergies.

## Features

- DOP-free
- No latex proteins, good alternative for those suffering from Type I allergies
- Suitable for food handling, except for fatty food
- Practical barrier for wide range of applications
- Good fit, feel and performance
- AQL 1,5

## Quality standards

- Complies with EN 455



## Technical information

Type	Powder-free, unsterile
Material	Vinyl
Colour	Transparent
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	S, M and L 240 mm
Impermeability as per EN 455-1	AQL 1,5
Durability, in original package if stored as per din 7716, ISO 2230	Min. 5 years

- Pcs/dispenser: 100
- Pcs/carton: 1000

**REF 210283S** S  
**REF 210283M** M  
**REF 210283L** L

# SELEFA® latex gloves

## Latex, smooth

Latex powdered examination gloves offer comfort and protection, a good combination of strength and elasticity. They perform well in heavy-duty situations and fit comfortably like a second skin. The only downside to natural rubber latex examination gloves is that being a protein based compound they may cause allergic reactions and irritation in some people – including anaphylactic reaction.

## Features

- Beaded cuff for easy donning and prevention of roll back
- Low level of extractable proteins and chemical residue
- Food safe, suitable for handling all foodstuffs
- AQL 1,5

## Quality standards

- Complies with EN 455



## Technical information

Type	Powdered, unsterile
Material	Latex
Colour	Natural white
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	XS, S, M, L and XL 240 mm
Impermeability as per EN 455-1	AQL 1,5
Durability, in original package if stored as per din 7716, ISO 2230	Min. 5 years

- Pcs/dispenser: 100
- Pcs/carton: 1000

**REF 210277XS** XS  
**REF 210277S** S  
**REF 210277M** M  
**REF 210277L** L  
**REF 210277XL** XL

# SELEFA® nitrile gloves

## Nitrile general use: blue and white

Comfortable, soft flexible powder-free nitrile gloves provide added safety in many applications. Designed with a special nitrile formulation they feel and fit like latex and allow full range of motion and excellent flexibility to minimise stress and fatigue. They do not contain natural rubber latex and are excellent alternative for those suffering from Type I allergies.

## Features

- 100% latex free – made from durable, protein and powder-free nitrile, eliminates the Type 1 allergic reaction associated with natural rubber latex
- Highly elastic and soft
- Rolled rim to facilitate easy donning
- Designed to give a natural rubber-like feel
- Textured in finger tips for enhanced grip in all situations
- Food safe, suitable for handling fatty foods incl. cooking, olive and cod liver oil
- Protects against a wide variety of chemicals
- AQL 1.5



## Quality standards

- Complies with EN 455:1–4 and EN 374:2

## Technical information

Type	Powder-free, unsterile
Material	Nitrile butadiene rubber (NBR)
Colour	Blue/white
Model	Ambidextrous, with rolled rim
Storage	Protect from heat, humidity, strong light and ozone
Size/overall length as per EN 455-2	XS, S, M, L and XL 240 mm
Impermeability as per EN 455-1	AQL 1,5
Durability, in original package if stored as per din 7716, ISO 2230	3 years

<b>Blue:</b>		<b>Blue:</b>	
<b>REF 210255XS</b>	XS	<b>REF 210253XS</b>	XS
<b>REF 210255S</b>	S	<b>REF 210253S</b>	S
<b>REF 210255M</b>	M	<b>REF 210253M</b>	M
<b>REF 210255L</b>	L	<b>REF 213253L</b>	L
<b>REF 210255XL</b>	XL	<b>REF 213253XL</b>	XL

Pcs/dispenser: 100  
Pcs/carton: 1000

Pcs/dispenser: 200  
Pcs/carton: 2000

<b>White:</b>		<b>White:</b>	
<b>210251XS</b>	XS	<b>210252XS</b>	XS
<b>210251S</b>	S	<b>210252S</b>	S
<b>210251M</b>	M	<b>210252M</b>	M
<b>210251L</b>	L	<b>210252L</b>	L
<b>210251XL</b>	XL	<b>210252XL</b>	XL

Pcs/dispenser: 100  
Pcs/carton: 1000

Pcs/dispenser: 200  
Pcs/carton: 2000

# SELEFA® wall dispensers

## Wall dispenser

Acrylic box dispenser for 1, 2 or 3 boxes of examination gloves (100 pcs/box or 150 pcs/box) or aprons.

- Pcs/unit pack: 1
- Pcs/transport carton: 1 x 10

**Mini**      **REF 215051**  
**Medi**      **REF 215052**  
**Maxi**      **REF 215053**

Suited for pcs of boxes

1 box

2 boxes

3 boxes



Cleaning: Mild soap solution and water or alternatively 70 % ethanol or similar.



## Shoe cover

### REF 211541

- Unsterile, PE film, blue
- Size: 41 x 15 cm
- Pcs/unit pack: 10
- Pcs/box: 100
- Pcs/carton: 2000



## Shoe cover

### REF 211545

- Unsterile, PE film, blue
- Size: 45 x 15 cm
- Pcs/unit pack: 10
- Pcs/box: 100
- Pcs/carton: 2000



## Shoe cover

### REF 211540

- Unsterile, PE film, white, heavy-duty
- Size: 40 x 15 cm
- Pcs/box: 50
- Pcs/carton: 1000



## Arm cover, elasticated cuff

### REF 210750

- Unsterile, PE film, white
- Size: 42 x 22 cm
- Pcs/unit pack: 10
- Pcs/box: 100
- Pcs/carton: 2000



## Arm cover, knitted cuff

### REF 210745

- Unsterile, PP nonwoven/PE film, nylon cuff, yellow
- Size: 21 x 43 cm
- Pcs/unit pack: 2
- Pcs/box: 50
- Pcs/carton: 400



# SELEFA® Protective apparel

## Visitor coat

### REF 210800

- Unsterile, PE film, white
- Size: 125 cm
- Pcs/unit pack: 1
- Pcs/box: 12
- Pcs/carton: 600



## Lab coat, velcro closure, 2 pockets, without cuff

### REF 210801

- Unsterile, PP nonwoven, white
- Size: 118 cm
- Pcs/unit pack: 1
- Pcs/carton: 50



## Beard cover

### REF 210855

- Unsterile, PP nonwoven, white
- Size: 48 x 24 cm
- Pcs/box: 100
- Pcs/carton: 1000



## Face mask 2-ply, paper

### REF 210451

- Unsterile, white
- Pcs/box: 100
- Pcs/carton: 10000



# SELEFA® Protective apparel

## Peaked cap

- Unsterile, PP nonwoven, cardboard peak
- White/white/blue peak
- Pcs/box: 100
- Pcs/carton: 1000

**REF 210763** M (50 cm)  
**REF 210767** L (55 cm)



## Peaked cap

- Unsterile, PP nonwoven, cardboard peak
- White/white/red peak
- Pcs/box: 100
- Pcs/carton: 1000

**REF 210764** M (50 cm)  
**REF 210768** L (55 cm)



## Snood peak with hairnet

**REF 210760**

- Unsterile, PP nonwoven
- White
- Size: XL
- Pcs/box: 100
- Pcs/carton: 1000



## Overall with hood

- Unsterile, spunbond PP nonwoven
- White
- Pcs/unit pack: 1
- Pcs/carton: 50

**REF 210810** M  
**REF 210811** L  
**REF 210812** XL



## Bed protection

Spreading of infection is prevented by protecting the patient mattress with a plastic mattress cover. Protection is especially important if the patient is incontinent or has exuding wounds. It is also good to protect the pillow from sweat and other secretions with a PE pillow cover. Isolation of patient bed is completed with draw sheet. OneMed assortment contains a large selection of draw sheets for different needs, also a polyethylene version.




# SELEFA® Bed protection

## Draw sheet Special, threadreinforced

### REF 500334

- White, unsterile
- Tissue/PE laminate, fluid barriers at long edges, reinforcement threads polyester
- Size: 80 x 140 cm
- Pcs/box: 25
- Pcs/carton: 100



## Draw sheet

### REF 500336

- White/blue, unsterile
- Tissue/PE laminate, fluid barriers at long edges
- Size: 80 x 140 cm
- Pcs/carton: 100



## Draw sheet Super, threadreinforced

### REF 500337

- White, unsterile
- Tissue/tissue/PE laminate, fluid barriers at long edges, reinforcement threads polyester
- Size: 80 x 140 cm
- Pcs/carton: 100



## Bed protection, roll

### REF 30680

- Unsterile, PE film
- Size: 90 x 130 cm
- Pcs/box: 50
- Pcs/carton: 500



## Mattress cover

### REF 030660

- Blue, unsterile, CPE film
- Size: 210 x 90 x 2 cm
- Pcs/unit pack: 1
- Pcs/box: 10
- Pcs/carton: 100



## Mattress cover

### REF 030655

- Blue, unsterile, PE film
- Size: 100 x 240 cm
- Pcs/box: 20
- Pcs/carton: 200



## Mattress cover

### REF 30670

- Transparent, unsterile, PE film
- Size: 75 x 160 cm
- Pcs/carton: 100

## Mattress cover

### REF 30665

- Transparent, unsterile, PE film
- Size: 95 x 230 cm
- Pcs/carton: 100

## Strether sheet

### REF 30780

- White, unsterile, viscose/polyester/cellulose nonwoven
- Size: 70 x 200 cm
- Pcs/carton: 100



## Examination table crepe

### REF 500510

- Natural white, unsterile
- Size: 50 x 200 cm
- Pcs/carton: 2 rolls



## Pillow case

### REF 30750

- White, unsterile, viscose/polyester/cellulose nonwoven
- Size: 50 x 70 cm
- Pcs/carton: 100



## Pillow case

### REF 030760

- Rose, unsterile, PE film
- Size: 60 x 77 cm
- Pcs/box: 100
- Pcs/carton: 1000



## Pillow case

### REF 030740

- Blue, unsterile, PE film
- Size: 40 x 50 cm
- Pcs/box: 100
- Pcs/carton: 1000



## Pillow bag, long side open

### REF 30751

- Transparent unsterile, PE film
- Size: 50 x 70 cm
- Pcs/box: 50
- Pcs/carton: 1000



## Pillow bag, short side open

### REF 30752

- Light green, unsterile, PE film
- Size: 50 x 75 cm
- Pcs/carton: 100
- Pcs/carton: 1500



## Bed sheet

### REF 30755

- Unsterile
- White, viscose/polyester/cellulose nonwoven
- Size: 140 x 240 cm
- Pcs/carton: 50



## Lifting sheet

### REF 30790

- Unsterile
- White PP nonwoven + transparent PE cast film
- Size: 140 x 240 cm
- Pcs/carton: 25



## Sampling & examining

The 2-layer injection swab is very absorbent. It is available both loose packed in box of 400 pcs and perforated in a roll. Each transport carton of roll-packed injection swabs contains one easy-to-use cardboard dispenser. For quick skin cleansing the swab assortment contains also an isopropyl alcohol skin cleansing swab.

The birch tree wooden spatulas – tongue depressor and cervical spatula – are packed into handy dispensers.

Two- or three-layer protective drapes are suitable for many uses. They easily provide a clean surface for sampling or for other procedures requiring clean and isolating surface.

The sampling & examining products also contains ear funnels (fitting on Heine otoscope) and a wide assortment of applicators.





## Protective drape

### REF 500505

- White/blue, unsterile
- Tissue/PE film laminate
- Size: 37 x 40 cm
- Pcs/box: 200
- Pcs/carton: 1200



## Protective drape

### REF 500506

- White/blue, unsterile
- Tissue/PE film laminate
- Size: 37 x 40 cm
- Pcs/box: 250
- Pcs/carton: 1500



## Protective drape

### REF 500507

- white, unsterile
- Tissue/PE film/tissue laminate
- Size: 37 x 40 cm
- Pcs/box: 200
- Pcs/carton: 1200



## Swab, perforated

### REF 22271

- White, unsterile
- Soft drylaid paper
- Size: 8 x 10 cm
- Pcs/carton: 500
- Pcs/carton: 6000



## Injection swab, perforated

### REF 22275

- White, unsterile
- Soft drylaid paper
- Size: 4 x 5 cm
- Pcs/carton: 1000
- Pcs/carton: 12000



## Injection swab, loose

### REF 22276

- White, unsterile
- Soft drylaid paper
- Size: 4 x 5 cm
- Pcs/box: 400
- Pcs/carton: 7200



## Injection swab, isopropyl alcohol 70 %

### REF 222680

- White, unsterile
- PP nonwoven, isopropyl alcohol
- Size: 3 x 3 cm
- Pcs/box: 200
- Pcs/carton: 4000



## Ear funnel

### REF 1201024

- Black, unsterile
- Size: 4 mm
- Pcs/box: 240
- Pcs/carton: 240



## Ear funnel

### REF 1201152

- Grey, unsterile
- Size: 2,5 mm
- Pcs/box: 200
- Pcs/carton: 200



## Ear funnel

### REF 1201022

- Black, unsterile
- Size: 2,5 mm
- Pcs/box: 240
- Pcs/carton: 240



## Ear funnel

### REF 1201154

- Grey, unsterile
- Size: 2,5 mm
- Pcs/box: 200
- Pcs/carton: 200



## Forceps

### REF 210500L

- White, unsterile
- Pcs/box: 50
- Pcs/carton: 300




# SELEFA® Sampling & examining

## Cotton tipped applicator, single tip

### REF 120783

- Unsterile
- Size: 15 cm
- Pcs/unit pack: 100
- Pcs/box: 2400
- Pcs/carton: 24000



## Cotton tipped applicator, double tip

### REF 120795

- Unsterile
- Size: 8 cm
- Pcs/unit pack: 100
- Pcs/box: 2400
- Pcs/carton: 24000



## Cotton tipped applicator

### REF 120788

- Sterile
- Size: 15 cm
- Pcs/unit pack: 2
- Pcs/box: 200
- Pcs/carton: 2000



## Cotton tipped applicator

### REF 120790

- Sterile
- Size: 15 cm
- Pcs/unit pack: 10
- Pcs/box: 250
- Pcs/carton: 2500



## Applicator sticks

### REF 120760

- Unsterile
- Size: 15 cm
- Pcs/box: 1000
- Pcs/carton: 50000



## Tongue depressor

### REF 140680

- Unsterile
- Size: 15 cm
- Pcs/box: 100
- Pcs/carton: 5000



## Cervical spatula

### REF 160298

- Unsterile
- Size: 18 cm
- Pcs/box: 100
- Pcs/carton: 5000



## Stethoscope

### REF 050947

- Unsterile
- Size: 60 cm
- Pcs/unit pack: 1
- Pcs/carton: 50



## Ear tips for stethoscope

### REF 509451

- Unsterile
- Size: 15 x 12 mm
- Pcs/unit pack: 2
- Pcs/carton: 1000



## Absorbent cotton wool in roll

### REF 221102

- Unsterile
- Size: 1 kg
- Pcs/unit pack: 1
- Pcs/carton: 16 rolls



## Towel

### REF 21201

- White, unsterile, drylaid paper
- Size: 40 x 40 cm
- Pcs/box: 50
- Pcs/carton: 300



## Patient hygiene & care

Traditional washing mitten is a basic product for patient's personal hygiene. SELEFA® range comprises of three different gloves: a drylaid paper version, drylaid paper/plastic film version and a thick needle punched nonwoven version. All are safe to use and comfortable for the patient. Of the same drylaid paper are made the patient wipes that can be used also for other cleaning and wiping. A thick drylaid paper towel is a good accessory for both washing and drying bed patients and for wiping and cleaning working tables and other surfaces.

Bathtub cover that is made of PE-plastic makes the bathing trolley new and personal for each patient at each bathing time.

Impermeable tie-on and tape fixation SELEFA® bibs protect the patient clothing well during meal times. The bibs are equipped with a pocket and are available in 2 sizes. In addition to traditional 2-layer model there is a 3-layer model. The materials of the soft bibs are crepe and PE-film.

Single packed toothbrushes with or without paste, nail and foot files and shower caps complete Patient hygiene & care assortment.



## Bib

### REF 500705

- Unsterile, laminate tissue/PE film
- Size: 37 x 60 cm
- Pcs/box: 100
- Pcs/carton: 900



## Bib

### REF 500706

- Unsterile, laminate tissue/PE film
- Size: 37 x 70 cm
- Pcs/box: 100
- Pcs/carton: 900



## Bib, triplex

### REF 500707

- Unsterile, laminate tissue/PE film/tissue
- Size: 37 x 70 cm
- Pcs/box: 100
- Pcs/carton: 600



## Bib, tape fixation

### REF 500708

- Unsterile, laminate tissue/PE film
- Size: 37 x 60 cm
- Pcs/box: 100
- Pcs/carton: 900

## Washglove

### REF 21101

- Unsterile, airlaid paper
- White
- Size: 16 x 22 cm
- Pcs/box: 100
- Pcs/carton: 1000



## Washglove with plastic lining

### REF 21102

- Unsterile, airlaid paper/PE film
- White
- Size: 16 x 22 cm
- Pcs/box: 100
- Pcs/carton: 1000



## Washglove extra soft

### REF 21103

- Unsterile, viscose/polyester nonwoven
- White
- Size: 16 x 22 cm
- Pcs/box: 100
- Pcs/carton: 1000



## Patient wipe

### REF 21104

- Unsterile, drylaid paper
- White
- Size: 20 x 20 cm
- Pcs/box: 100
- Pcs/carton: 1400





## Napkin

### REF 210101

- Unsterile, tissue
- Size: 33 x 33 cm
- Pcs/box: 100
- Pcs/carton: 4600



## Cleaning wipe

### REF 21202

- Unsterile, nonwoven, viscose/polyester/cotton/cellulose fibres
- Blue
- Size: 50 x 60 cm
- Pcs/carton: 400



## Toothbrush with paste

### REF 210645

- Unsterile
- Pcs/pack: 1
- Pcs/box: 100
- Pcs/carton: 1000



## Toothbrush without paste

### REF 210640

- Unsterile
- Pcs/pack: 1
- Pcs/box: 100
- Pcs/carton: 1000



## Foot file

### REF 210628

- Unsterile
- Pcs/box: 100
- Pcs/carton: 2500



## Nail file

### REF 210625

- Unsterile
- Pcs/box: 100
- Pcs/carton: 5000



## Shower cap

### REF 030635

- Unsterile, PE film
- Pcs/box: 100
- Pcs/carton: 2000

## Flexible drinking straw

### REF 030344

- Unsterile
- Size: 21 cm
- Pcs/unit pack: 1
- Pcs/box: 200
- Pcs/carton: 5000



## Triangular sling

### REF 222415

- Unsterile, white
- Size: 95 x 95 x 135 cm
- Pcs/unit pack: 1
- Pcs/carton: 50



## Triangular sling

### REF 222410

- Unsterile, natural white
- Size: 90 x 90 x 130 cm
- Pcs/carton: 600

## Basin bag

### REF 210620S

- Unsterile, transparent, PE film
- Size: 65 x 40 cm
- Pcs/box: 250
- Pcs/carton: 5000



## Basin bag

### REF 210621

- Unsterile, transparent, PE film
- Size: 70 x 47 cm
- Pcs/box: 50
- Pcs/carton: 1000

## Bath tub cover, roll

### REF 30701

- Unsterile, PE film
- Transparent
- 135 m x 135/270 cm
- Pcs/carton: 1 roll

## Wound care

SELEFA® nonwoven swabs and sponges are cost-effective alternatives to traditional cotton gauze. They feature superior absorption capacity, softness and low fibre release. The raw material is high quality viscose/polyester nonwoven.

SELEFA® gauze swabs, sponges and tamponades are available in large number of sizes and packaging configurations. They are made of 17-thread gauze that complies with the stipulations of the European Pharmacopoeia.

SELEFA® products are packed into easy-to-use packaging.



# SELEFA® Wound care

REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
-----	---------------------	------	---------------	---------------	------------

## Gauze swabs, unsterile, 8-ply packed in paper bag

<b>22219513</b>	Gauze swabs	5 x 5 cm	-	100	1000
<b>22219524</b>	Gauze swabs	7,5 x 7,5 cm	-	100	1000
<b>22219535</b>	Gauze swabs	10 x 10 cm	-	100	1000
<b>22219546</b>	Gauze swabs	10 x 20 cm	-	100	2000

22219513



## Gauze swabs, unsterile, 12-ply packed in paper bag

<b>22219512</b>	Gauze swabs	5 x 5 cm	-	100	1000
<b>22219523</b>	Gauze swabs	7,5 x 7,5 cm	-	100	1000
<b>22219534</b>	Gauze swabs	10 x 10 cm	-	100	1000
<b>22219545</b>	Gauze swabs	10 x 20 cm	-	100	1000

## Gauze swabs, presterilized, unsterile, 8-ply, packed in dispenser

<b>22219511</b>	Gauze swabs	5 x 5 cm	-	120	10800
<b>22219522</b>	Gauze swabs	7,5 x 7,5 cm	-	120	7200
<b>22219533</b>	Gauze swabs	10 x 10 cm	-	120	7200

22219533



## Gauze swabs, presterilized, unsterile, 12-ply, packed in dispenser

<b>22219544S</b>	Gauze swabs	10 x 20 cm	-	50	1000
------------------	-------------	------------	---	----	------

## Gauze swabs, presterilized, unsterile, 16-ply, packed in dispenser

<b>22219555</b>	Gauze swabs	10 x 10 cm	-	60	3600
-----------------	-------------	------------	---	----	------

222250T



280500

## Triangular and flat gauze sponges, presterilized, unsterile

<b>222250T</b>	Triangular gauze sponges	30 x 30 x 30 mm	100	-	10000
<b>280500</b>	Flat gauze sponges	8 x 8 mm	100	1000	10000
<b>222260S</b>	Flat gauze sponges	30 x 30 mm	100	100	10000

222260S



## Gauze squares, presterilized, unsterile

<b>22217815</b>	Gauze squares	15 x 15 cm	-	500	10000
<b>22217820</b>	Gauze squares	20 x 20 cm	-	500	10000
<b>22217830</b>	Gauze squares	30 x 30 cm	-	500	5000
<b>22217840</b>	Gauze squares	40 x 40 cm	-	500	4000
<b>22217850</b>	Gauze squares	50 x 50 cm	-	250	2500

22217820



# SELEFA® Wound care

REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
-----	---------------------	------	---------------	---------------	------------

## Gauze round sponges, unsterile

<b>222253</b>	Gauze round sponges	Ø 30 mm	-	100	2000
<b>222256</b>	Gauze round sponges	Ø 70 mm	-	100	1000

## Gauze round sponges, presterilized, unsterile

<b>222250S</b>	Gauze round sponges	Ø 30 mm	-	200	3600
<b>222250M</b>	Gauze round sponges	Ø 40 mm	-	150	2700
<b>222250L</b>	Gauze round sponges	Ø 50 mm	-	100	1800
<b>222250XL</b>	Gauze round sponges	Ø 65 mm	-	100	1400

222250L



2222722



2222745



## Gauze swabs, sterile, 8-ply, 2 or 5 pcs/bag

<b>2222702</b>	Gauze swabs	5 x 5 cm	2	120	2400
<b>2222722</b>	Gauze swabs	7,5 x 7,5 cm	2	120	2400
<b>2222742</b>	Gauze swabs	10 x 10 cm	2	120	2160
<b>2222705</b>	Gauze swabs	5 x 5 cm	5	125	2500
<b>2222725</b>	Gauze swabs	7,5 x 7,5 cm	5	125	2500
<b>2222745</b>	Gauze swabs	10 x 10 cm	5	125	2250



## Gauze swabs, sterile, 12-ply, 2 or 5 pcs/bag

<b>222291</b>	Gauze swabs	5 x 5 cm	2	80	1440
<b>222290</b>	Gauze swabs	7,5 x 7,5 cm	2	80	1600
<b>222289</b>	Gauze swabs	10 x 10 cm	2	80	1280
<b>2222762</b>	Gauze swabs	10 x 20 cm	2	40	720
<b>222296</b>	Gauze swabs	5 x 5 cm	5	80	1600
<b>222295</b>	Gauze swabs	7,5 x 7,5 cm	5	80	1600
<b>222294</b>	Gauze swabs	10 x 10 cm	5	80	1280
<b>2222765</b>	Gauze swabs	10 x 20 cm	5	40	720

222295



# SELEFA® Wound care

REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
-----	---------------------	------	---------------	---------------	------------

## Gauze swabs, sterile, 8-ply, special sizes

<b>222298</b>	Gauze swabs	20 x 80 cm	5	25	150
<b>222299</b>	Gauze swabs	10 x 60 cm	5	65	390
<b>222300</b>	Gauze swabs	10 x 10 cm	25	150	2400

## Gauze, round, triangular and flat sponges, sterile

<b>22225030</b>	Round gauze sponges	Ø 30 mm	5	75	1500
<b>22225040</b>	Round gauze sponges	Ø 40 mm	5	75	1350
<b>22225050</b>	Round gauze sponges	Ø 50 mm	5	50	900
<b>222254</b>	Round gauze sponges	Ø 70 mm	5	100	800
<b>222250TS</b>	Triangular gauze sponges	30 x 30 x 30 mm	5	250	5000
<b>22226030</b>	Flat gauze sponges	30 x 30 mm	5	250	5000



## Nonwoven swabs, unsterile, 4-ply, packed in paper bag

<b>22319512</b>	Nonwoven swabs	5 x 5 cm	-	100	6000
<b>22319523</b>	Nonwoven swabs	7,5 x 7,5 cm	-	100	7500
<b>22319534</b>	Nonwoven swabs	10 x 10 cm	-	100	5400
<b>22319545</b>	Nonwoven swabs	10 x 20 cm	-	100	2700



## Nonwoven swabs, presterilized, unsterile, 4-ply, packed in box

<b>22319511</b>	Nonwoven swabs	5 x 5 cm	-	150	13500
<b>22319522</b>	Nonwoven swabs	7,5 x 7,5 cm	-	150	9000
<b>22319533</b>	Nonwoven swabs	10 x 10 cm	-	150	9000
<b>22319544</b>	Nonwoven swabs	10 x 20 cm	-	75	1500



# SELEFA® Wound care

REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
<b>Nonwoven round sponges, presterilized, unsterile</b>					
<b>223250S</b>	Nonwoven round sponges	Ø 30 mm	-	200	3600
<b>223250M</b>	Nonwoven round sponges	Ø 40 mm	-	150	2700
<b>223250L</b>	Nonwoven round sponges	Ø 50 mm	-	100	1800
<b>223250T</b>	Triangular nonwoven sponges	30 x 30 x 30 mm	100	1000	10000



## Nonwoven swabs, sterile, 4-ply, packed 2 or 5 pcs/bag

<b>2232702</b>	Nonwoven swabs	5 x 5 cm	2	140	2800
<b>2232722</b>	Nonwoven swabs	7,5 x 7,5 cm	2	140	2800
<b>2232742</b>	Nonwoven swabs	10 x 10 cm	2	140	2520
<b>2232762</b>	Nonwoven swabs	10 x 20 cm	2	70	1260
<b>2232705</b>	Nonwoven swabs	5 x 5 cm	5	150	3000
<b>2232725</b>	Nonwoven swabs	7,5 x 7,5 cm	5	150	3000
<b>2232745</b>	Nonwoven swabs	10 x 10 cm	5	150	2700
<b>2232765</b>	Nonwoven swabs	10 x 20 cm	5	40	720



## Nonwoven round sponges, sterile, packed 5 pcs/bag

<b>22325030</b>	Round nonwoven sponges	Ø 30 mm	5	75	1500
<b>22325040</b>	Round nonwoven sponges	Ø 40 mm	5	75	1350
<b>22325050</b>	Round nonwoven sponges	Ø 50 mm	5	50	900
<b>223259</b>	Round nonwoven sponges	Ø 70 mm	5	50	400
<b>223250TS</b>	Triangular nonwoven sponges	30 x 30 x 30 mm	5	250	5000



## Nonwoven round sponges, unsterile

<b>223253</b>	Nonwoven round sponges	Ø 30 mm	-	100	2000
<b>223256</b>	Nonwoven round sponges	Ø 70 mm	-	100	1000

# SELEFA® Wound care

REF	Product description	Size	Pcs/unit pack	Pcs/dispenser	Pcs/carton
<b>Tamponing bandages, sterile</b>					
222463	Tamponing bandage	1 x 50 cm	1	10	120
222464	Tamponing bandage	1 x 100 cm	1	10	120
222465	Tamponing bandage	1 x 200 cm	1	10	120
222466	Tamponing bandage	1 x 500 cm	1	10	120
222467	Tamponing bandage	2 x 50 cm	1	10	120
222468	Tamponing bandage	2 x 100 cm	1	10	120
222469	Tamponing bandage	2 x 200 cm	1	10	120
222470	Tamponing bandage	2 x 500 cm	1	10	120
222471	Tamponing bandage	5 x 50 cm	1	10	120
222472	Tamponing bandage	5 x 100 cm	1	10	120
222473	Tamponing bandage	5 x 200 cm	1	10	120
222474	Tamponing bandage	5 x 300 cm	1	10	120
222475	Tamponing bandage	5 x 500 cm	1	10	120
222476	Tamponing bandage	10 x 100 cm	1	10	120
222477	Tamponing bandage	10 x 300 cm	1	10	120



222465



222473



## Tamponing bandage, x-ray detectable, sterile

226055SX	Tamponing bandage	Ø 55 mm	1	25	100
----------	-------------------	---------	---	----	-----

226055SX






# SELEFA® Surgical tapes

## Surgical tape, nonwoven with holder

A latex-free, breathable nonwoven tape.  
Gentle to the skin yet adheres well.  
An economical, general purpose surgical tape.

### REF 2405021

- White
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720



### REF 2405022

- White
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360



### REF 2405023

- White
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180



## Surgical tape, nonwoven

A latex-free, breathable nonwoven tape.  
Gentle to the skin yet adheres well.  
An economical, general purpose surgical tape.

### REF 2402021

- White
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720

### REF 2402022

- White
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360



### REF 2402023

- White
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180

## Surgical tape, Silk

Latex-free, white, silk like high strength surgical tape.  
Bi-directional tearing for easy tailoring.

### REF 2402011

- White
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720

### REF 2402012

- White
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360



### REF 2402013

- White
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180

## Surgical tape, PE-film

A transparent medical pe-tape.  
Perforated, that tears easily.

### REF 2402031

- Transparent
- Size: 1,25 cm x 9,1 m
- Pcs/dispenser: 24
- Pcs/carton: 720

### REF 2402032

- Transparent
- Size: 2,5 cm x 9,1 m
- Pcs/dispenser: 12
- Pcs/carton: 360



### REF 2402033

- Transparent
- Size: 5 cm x 9,1 m
- Pcs/dispenser: 6
- Pcs/carton: 180

## Injection plasters

### Adhesive Plaster, Nonwoven

#### REF 240100F

- Skin
- Size: 24 x 70 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

#### REF 240101F

- Skin
- Size: 10 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

### Adhesive Plaster, Plastic

#### REF 240100FT

- Transparent
- Size: 24 x 70 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

#### REF 240101FT

- Transparent
- Size: 10 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



## Round plasters

### Adhesive Plaster, Nonwoven

#### REF 240102F

- Skin
- Size: 38 x 38 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



#### REF 240103F

- Skin
- Size: 20 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



### Adhesive Plaster, Plastic

#### REF 240102FT

- Transparent
- Size: 38 x 38 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200



#### REF 240103FT

- Transparent
- Size: 20 x 40 mm
- Pcs/dispenser: 200
- Pcs/carton: 10200

## Childrens plasters

### Adhesive Plaster, Plastic "Pinocchio"

#### REF PHA20PML

- Motive
- Size: 20 x 70 mm
- Pcs/dispenser: 240



